

## ENDOCUFF VISION™ Compatibility Chart\*

OLYMPUS		
Series	Scope Model	ENDOCUFF VISION™ Model
EVIS EXERA III	CF-HQ190L/I	ARF120 (Green)
	CF-H190L/I	ARF120 (Green)
	PCF-H190DL/I	ARF110 (Blue)
	PCF-H190L/I	ARF130 (Purple)
	PCF-PH190L/I	Not compatible
	CF-H185L/I	ARF110 (Blue)
EVIS EXERA II	CF-H180AL/I	ARF120 (Green)
	CF-H180DL/I	ARF140 (Orange)
	CF-Q180AL/I	ARF120 (Green)
	PCF-H180AL/I	ARV130 (Purple)
	PCF-Q180AL/I	ARV130 (Purple)
	CF-Q165L/I	ARF120 (Green)
EVIS EXERA	CF-Q160AL/I	ARF120 (Green)
	CF-Q160DL/I	ARF120 (Green)
	CF-Q160L/I	ARF120 (Green)
	CF-Q160S	ARF120 (Green)
	PCF-160AL/I	ARV130 (Purple)
	CF-Q145L/I	ARF110 (Blue)
EVIS 140	CF-Q140L/I	ARF120 (Green)
	CF-140L/I	ARF120 (Green)
	CF-140S	ARF120 (Green)
	CF-1T140L/I	ARF120 (Green)
	PCF-140L/I	Not compatible
EVIS LUCERA ELITE	CF-HQ290L/I	ARF120 (Green)
	CF-H290L/I	ARF120 (Blue)
	PCF-H290L/I	ARV130 (Purple)
EVIS LUCERA SPECTRUM	CF-H260AL/I	ARF120 (Green)
	CF-H260DL/I	ARF120 (Green)
	CF-H260AZL/I	ARF120 (Green)
	CF-Q260AL/I	ARF110 (Blue)
	CF-Q260DL/I	ARF110 (Blue)
	PCF-Q260AL/I	ARV130 (Purple)
	PCF-PQ260L/I	Not compatible
	PCF-Q260AZI	ARV130 (Purple)
EVIS LUCERA	CF-Q240AL/I	ARF120 (Green)
	CF-Q240L/I	ARF120 (Green)
	CF-Q240ZL/I	ARF120 (Green)
	CF-240AL/I	ARF110 (Blue)
	CF-240DL/I	ARF110 (Blue)
	CF-240L/I	ARF110 (Blue)
OPTERA	CF-H170L/I	ARF120 (Green)

\*list is limited to devices whose compatibility to ENDOCUFF VISION (TM) has been tested by manufacturer of ENDOCUFF VISION Boddingtons during latest release June 2017

**FUJI**

Scope Model	ENDOCUFF VISION™ Model
EC-760R-V/L	ARV110 (Blue)
EC-760ZP-VM/VL	ARV110 (Blue)
EC-600WM/WI/WL	ARV110 (Blue)
EC-590WM4/WI4/WL4/WI/WLY	ARV110 (Blue)
EC-590ZW3 M/L	ARV120 (Green)
EC-530FI/FL/FM	ARV120 (Green)
EC-530WM3/WI3/WL3/WI	ARV120 (Green)
EC-450 Adult Series	ARV120 (Green)
EC-200MR	ARV120 (Green)
EC-250 Adult Series	ARV120 (Green)
EC-201WI	ARV120 (Green)
EC-530 MT/IT/LT	ARV120 (Green)
EC-530 DM/DL	ARV120 (Green)
EC-600 ZW/M/L	ARV120 (Green)
EC-530MP/LP	ARV130 (Purple)
EC-530LS	ARV130 (Purple)
EC-450 Paediatric Series	ARV130 (Purple)
EC-450LP/LS5	ARV130 (Purple)
EC-250 Intermediate Series	ARV130 (Purple)

**PENTAX**

EC38-i10M/M2/F/F2/L/NF	ARV110 (Blue)
EC-3890TMK/MK2/FK/FK2/LK	ARV110 (Blue)
EC-3890MZI/FZI/LZI/SI2/TFK	ARV110 (Blue)
EC-3890FI/LI/Fi2(HD)/Mi2	ARV110 (Blue)
EC-3890TLK	ARV110 (Blue)
EC-3970LK	ARV120 (Green)
EC-3880FK2	ARV120 (Green)
EC-3870LK/LZK/TLK/FK2/MK2	ARV120 (Green)
EC-3840FK2	ARV120 (Green)
EC-3830LK/P/FK2	ARV120 (Green)
EC-380LKP/FKP/FK2P/MK2P	ARV120 (Green)
EC-38LV	ARV120 (Green)
EC-3940FK	ARV130 (Purple)
EC-3490 FI	ARV130 (Purple)
EC-3490 MK/FK/LK/LI/TLI/TFI	ARV130 (Purple)
EC-3470LK	ARV130 (Purple)
EC-34-i10M/F/L	ARV130 (Purple)
EC-3490TLF/MI	Not compatible
EC-2990LI	Not compatible

**STORZ**

13924-PKS	ARV120 (Green)
13925-PKS	ARV120 (Green)
13907-PKS	ARV130 (Purple)
13906-PKS	ARV130 (Purple)
13905-PKS	ARV140 (Orange)
13904-PKS	ARV140 (Orange)

\*list is limited to devices whose compatibility to ENDOCUFF VISION (TM) has been tested by manufacturer of ENDOCUFF VISION Boddingtons during latest release June 2017

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.



Distributed by

**OLYMPUS EUROPA SE & CO. KG**  
 Postbox 10 49 08, 20034 Hamburg, Germany  
 Wendenstrasse 14-18, 20097 Hamburg, Germany  
 Phone: +49 40 23773-0, Fax: +49 40 233765  
 www.olympus-europa.com

**ENDOCUFF VISION™**  
 Trademark holder  
 Arc Medical Design Ltd  
 43 Park Place  
 Leeds LS1 2RY  
 United Kingdom